

Third Millennium Laser Therapy.

MLS[®]
Multiwave
Locked System




ASV
A S A L A S E R
Research and Therapeutic Solutions

What makes our laser different from the rest? Three letters: MLS®.

MLS® (Multiwave Locked System) combines and synchronises specific continuous and pulsed emissions with different wavelengths. This synchronisation is produced thanks to a patented control system, which allows modulation of emission type ac-

cording to requirements. The synchronised and controlled laser energy reaches deep into the tissue to create a more powerful therapeutic anti-inflammatory, anti-edema and analgesic effect speeding up the healing process at the same time.

- ▶ Synergy and strengthening of the anti-inflammatory, analgesic, anti-edema effects
- ▶ Fast pain reduction
- ▶ Improvement of local microcirculation
- ▶ Short treatment times and few treatments required (usually from 3 to 10 sessions)



MLS® has been scientifically proven on the ASAcampus.

The efficiency and safety of MLS® Laser Therapy has been proved scientifically on the ASAcampus, the Joint Laboratory at the Department of Biomedical and Clinical Experimental Science at the University of Florence - Italy - directed by Dr. Monica Monici.

LEARN MORE ABOUT ASACAMPUS



▶ Mphi



▶ Mphi 75



▶ Mphi 5



▶ Mphi 75.5



▶ M6



▶ M8





CORPORATE HEADQUARTERS / REGISTERED OFFICE

Via Galileo Galilei, 23 / 36057 Arcugnano (VI) - Italy

T +39 0444 28 92 00 / F +39 0444 28 90 80

asalaser@asalaser.com

RESEARCH DIVISION / BRANCH

Joint Laboratory Department of Experimental and Clinical

Biomedical Sciences University of Florence

Viale G. Pieraccini, 6 / 50139 Firenze - Italy

asacampus@asalaser.com

asalaser.com

